


# Loading Samples into the XPS Endstation

## Component Reference:

- **A:** Nitrogen gas valve
- **B:** Sample rack knob
- **C:** Gate valve between main chamber and load-lock
- **D:** Roughing pump valve knob
- **E:** Turbo pump power cord plug



# Part 1: Loading Sample Plates into the Load-Lock

## 1. Ensure the Gate Valve is Closed

Before venting the load-lock with nitrogen gas, **confirm that the gate valve (C) between the main chamber and load-lock is fully closed.**



## 2. Turn Off the Ion Gauge (IMG2)

a) On the ion gauge controller, use the drop-down arrow to select **IMG2**, which shows the **load-lock pressure**.

b) Press the “-” button until “off” appears, then press “OK” to turn it off. The IMG2 pressure reading should drop to **zero**.

(2-a)



(2-b)



### 3. Unplug the Turbo Pump

Disconnect the **LL-labelled power cord** of the turbo pump from the electric outlet.



### 4. Close the Roughing Pump Valve

Turn the **roughing pump valve (D)** clockwise to close it.





### 5. Open the Nitrogen Gas Cylinder

Turn the cylinder valve **counter-clockwise** to fully open it. Check the **Swagelok gauge** to ensure the pressure is around **20 psi**.



## 6. Introduce Nitrogen Gas

Slowly open the **nitrogen valve (A)** counter-clockwise (6-a). Use the **TCG controller** (6-b) to monitor the pressure and introduce nitrogen at a rate of **1–2 Torr per second**.



## 7. Vent the Load-Lock

Loosen the **load-lock door knob** (counter-clockwise) and wait for the TCG controller to show a pressure of approximately **700 Torr**.



## 8. Open the Load-Lock Door

Once pressure exceeds 700 Torr, use a screwdriver to **gently open the load-lock door**. **Keep nitrogen flowing** while the door is open.

## 9. Load the Sample Plates


Turn the **sample rack knob (B)** to extend the rack. Load up to **5 sample plates**, ensuring each is securely fixed.



**Note:** Before closing the load-lock door, **clean the door and O-ring** using a Kimwipe lightly moistened with **methanol**.

## 10. Close the Load-Lock

Fully close the **nitrogen valve (A)** (clockwise), then close and **finger-tighten** the load-lock door.



## 11. Restart the Pumping Process

- a) Close the **nitrogen gas cylinder** (clockwise).
- b) Open the **roughing pump valve (D)** (counter-clockwise).

(11-a)



(11-b)



### 12. Wait for Initial Vacuum

Monitor the TCG controller until pressure drops below **100 mTorr**.




### 13. Reconnect the Turbo Pump

Plug in the turbo pump's power cord.

- a) A **green light** will flash on the turbo pump, then stabilize.
- b) Wait until the TCG controller reads approximately **-0.0 mTorr**.

(13-a)



(13-b)



## 14. Turn On the Ion Gauge (IMG2)

- Select **IMG2** on the controller.
- Press the “+” button, then “OK” to turn it on.
- Wait until IMG2 shows a pressure of  **$10^{-8}$  Torr or lower** before proceeding.




---

## Part 2: Transferring a Sample Plate to the Main Chamber

### 1. Check Load-Lock Vacuum

Ensure **IMG2** pressure reads  $\leq 10^{-8}$  Torr.




### 2. Open the Gate Valve

Open the **gate valve (C)** (counter-clockwise) between the **main chamber and load-lock**.



### 3. Confirm Main Chamber Vacuum

Verify **IMG1** pressure remains around  $10^{-9}$  to  $10^{-10}$  Torr.



### 4. Use the Transfer Arm

Use the **transfer arm** to retrieve a sample plate. Insert it into the **main chamber** carefully, keeping the arm steady and avoiding rotation.



**Note:** Only **one sample plate** may be transferred at a time.

### 5. Mount the Sample

Align the sample holder using the **sample rack knob (B)** and gently insert the plate into the rack.



## **6. Retract and Close**

Return the transfer arm to its original position and **close the gate valve** between the main chamber and load-lock.

## **7. Ready for Measurement**

Your sample is now loaded, and the system is ready for measurement.

---

**Last Updated:** May 01, 2025